



WHITE PAPERS & APPLICATION NOTES



DOWNLOAD FREE WHITE PAPERS & APPLICATION NOTES

Highly Efficient Tunable Filter for Broadband Light Source

A new, simple and compact approach for tunable wavelength filtering is introduced with a device called Flexible Wavelength Selector. This new device combines the flexibility of a monochromator with the imaging functionality of a thin-film interference filter, and is ideal for use in both illumination and image acquisition in optical microscopes. The main features include high transmission efficiency, circular round aperture of up to 10 mm and independent tuning of center wavelength and bandwidth.

DOWNLOAD NOW

SLI
Light Done Right

Highly efficient tunable filter for broadband light source

Jun Hee Kang, Spectrolight Inc. 19800 MacArthur Blvd., Suite 300, Irvine, CA 92612

Abstract
A new simple and compact approach to tunable wavelength filtering combines the flexibility of a monochromator with the imaging functionality of a thin-film filters, and is ideal for use in both illumination and image acquisition.

Traditional tools for wavelength filtering, such as monochromators and filter wheels, all have some type of limitation that often compromises their use in scientific and industrial settings. However, a new type of wavelength filtering device called Flexible Wavelength Selector (FWS) now offers a combination of advantages for both illumination and image filtering. Specifically, FWS devices combine the wavelength tunability and bandwidth control of a monochromator with the circular uniform aperture of a filter wheel. These devices are highly efficient, simple, and cost-effective. In this article, we describe this innovative product and show its application with broadband lamp.

Figure 1. FWS Poly

Methods for filtering light according to wavelength can be divided into three groups: filters (cut-off, dichroic, etc.), monochromators, and newer "niche" technologies. Each method has its advantages.

Filters. In the detection part of a microscope, filters and dichroic beamsplitters are the commonly used means of wavelength filtering just prior to the camera or eyepiece. Similarly, these filters can be used with a broadband light source for wavelength-selective illumination instead of an LED or laser. The filter is typically based on thin film dielectrics sometimes incorporating colored glass for extra wavelength blocking. The limitation of filters and filter

19800 MacArthur Blvd., Suite 300, Irvine, CA 92612 | Tel: 949-450-7700
www.spectrolight.com

SPECTROLIGHT Inc.

Sponsored by



PHOTONICS MEDIA

Visit Photonics Media to download other white papers and learn more about the latest developments in lasers, imaging, optics, biophotonics, machine vision, spectroscopy, microscopy, photovoltaics and more.

www.photonics.com/WhitePapers.aspx

We respect your time and privacy. You are receiving this email because you are a Photonics Media subscriber, and/or a member of our website, Photonics.com. You may use the links below to manage your subscriptions or contact us.

Questions: info@photonics.com

[Unsubscribe](#) | [Subscribe](#) | [Subscriptions](#) | [Privacy Policy](#) | [Terms and Conditions of Use](#)

Photonics Media, 100 West St., PO Box 4949, Pittsfield, MA 01202-4949
© 1996 - 2019 Laurin Publishing. All rights reserved. Photonics.com is Registered with the U.S. Patent & Trademark Office.
Reproduction in whole or in part without permission is prohibited.