



WHITE PAPERS & APPLICATION NOTES



DOWNLOAD FREE WHITE PAPERS & APPLICATION NOTES

Study of Fluorescence Quenching Kinetics Using Stopped Flow

In this application note the kinetics of NATA quenching by QBS are determined using a Photoluminescence Spectrometer equipped with a stopped-flow accessory.

DOWNLOAD NOW

APPLICATION NOTE
Study of Fluorescence Quenching Kinetics Using Stopped-Flow
 AN_P31: 12 Dec. 2017, Stuart Thomson, Anna Gokarny

Introduction
 The understanding of reaction kinetics is an important area of study across biology and chemistry. Through the study of reaction kinetics the reaction rate, reaction order and the underlying molecular mechanisms can be determined. A powerful method for determining the kinetics of a reaction is by monitoring the concentration of the reactants or products over time using absorption or fluorescence spectroscopy. For slow reactions, the reactants can simply be mixed by hand or using a magnetic stirrer. However, for fast reactions where the lifetime of the reaction is comparable with the mixing time, more specialised techniques must be used such as stopped-flow.

Methods and Materials
 N-acetylcysteine (NAC) and N-bromosuccinimide were dissolved in phosphate buffered saline (PBS) at pH 7.3. Three NATA solutions were prepared with concentrations of 5 µM, 10 µM and 20 µM. The NBS quencher was kept in excess with concentrations of 20 µM, 100 µM and 200 µM respectively. The absorption spectrum of NATA was recorded using the transmission detector of the PFS Spectrofluorimeter with an excitation bandwidth of 3 nm. Fluorescence measurements were performed using the FLS1000 Photoluminescence Spectrometer equipped with double excitation and emission monochromators, photomultiplier tube detector (PMT-NBS), a 450 W Xe lamp and the N-MORMM stopped-flow accessory.

The emission spectrum of NATA was recorded at an excitation wavelength of 280 nm using excitation and emission bandwidths of 3 nm and 1.5 nm respectively. For fluorescence kinetic measurements the NATA and NBS solutions were loaded into the A & B syringes of the N-MORMM stopped-flow accessory. Syringe C was filled with PBS buffer to serve as a reference. The TTL output trigger of the FLS1000 was used to initiate the N-MORMM injection and the total stop volume was set to 400 µl. The kinetic spectra were recorded using an excitation wavelength of 280 nm with a 5 nm bandpass, an emission wavelength of 360 nm with a 5 nm bandpass, and a temporal resolution of 0.2 ns.

Results Discussion
 The absorption and emission spectra of NATA solution were recorded and shown in Figure 1. NATA has an absorption maximum at 280 nm and an emission maximum close to 360 nm and these were therefore chosen as the excitation and emission wavelengths for the subsequent fluorescence kinetic measurements.

A fluorescence kinetic of the quenching of 10 µM NATA by 100 µM NBS is shown in Figure 2 (red). The fluorescence quenching kinetic was analysed using the FLS1000 Fluorescence operating software. The data was fit well using a single monoexponential decay with a time constant of 28 ns. The fluorescence kinetic was recorded with a high temporal resolution.

Figure 1: Absorption (blue) and emission (red) spectra of NATA in PBS buffer.

Figure 2: Oxidation of N-acetylcysteine (NAC) by N-bromosuccinimide (NBS).

Sponsored by



More Application Notes from this Sponsor

- Electroluminescence and Photoluminescence Spectroscopy of a Phosphorescent Organic Light Emitting Diode (PhOLED)
- Photophysical Characterisation of Perovskite Quantum Dots
- Charge Carrier Recombination Dynamics of Semiconductor Photocatalyst

PHOTONICS MEDIA

Visit Photonics Media to download other white papers and learn more about the latest developments in lasers, imaging, optics, biophotonics, machine vision, spectroscopy, microscopy, photovoltaics and more.

www.photonics.com/WhitePapers.aspx

We respect your time and privacy. You are receiving this email because you are a Photonics Media subscriber, and/or a member of our website, Photonics.com. You may use the links below to manage your subscriptions or contact us.

Questions: info@photonics.com

[Unsubscribe](#) | [Subscribe](#) | [Subscriptions](#) | [Privacy Policy](#) | [Terms and Conditions of Use](#)

Photonics Media, 100 West St., PO Box 4949, Pittsfield, MA 01202-4949

© 1996 - 2019 Laurin Publishing. All rights reserved. Photonics.com is Registered with the U.S. Patent & Trademark Office. Reproduction in whole or in part without permission is prohibited.