



# WHITE PAPERS & APPLICATION NOTES




**DOWNLOAD FREE WHITE PAPERS & APPLICATION NOTES**

## Kinetics of Photocatalysis Reactions Studied by Transient Absorption Spectroscopy

Photocatalysis is the rate increase of a chemical reaction by light, often in the presence of a catalyst that starts the reaction upon irradiation. Photocatalysts are typically semiconducting metal oxides which are employed as particles in solution. Discover how Edinburgh Instruments LP980 Spectrometer was used in this application note.

**DOWNLOAD NOW**

**APPLICATION NOTE**  
**Kinetics of Photocatalysis Reactions Studied by Transient Absorption Spectroscopy**  
 ANI\_P90, Maria Tessa



**Introduction**  
 Photocatalysis is the rate increase of a chemical reaction by light, often in the presence of a catalyst that starts the reaction upon irradiation. Photocatalysts are typically semiconducting metal oxides such as  $ZnO$ ,  $Fe_2O_3$  or  $TiO_2$  which are employed as particles in solution. When absorbing light, these materials are able to generate electrons and holes which go on to react with chemical species on their surface.

The photocatalysis process is illustrated in Figure 1. If a photon of greater energy than the bandgap of the semiconductor is absorbed by the photocatalyst, it excites an electron (e) from the valence band into the conduction band. This process generates a hole (h) or positive charge in the valence band. These charge carriers are mobile, albeit short lived, and can excite through multiple pathways. The electron-hole pair can recombine emitting energy after a short time, but the charge carriers can also be trapped at specific sites in the material. If they do not recombine, electrons and holes can react with species on the surface of the photocatalyst. Valence band electrons can combine with electron-acceptors in solution, and conduction band holes may react with electron donors. This is often the first step in a series of charge-transfer reactions in solution that lead to a desired product.




Figure 1 Charge carrier generation in photocatalysis. Conduction band electrons react with electron-acceptors (EA) and valence band holes react with electron donors (ED).




Figure 2 Edinburgh Instruments LP980 Transient Absorption Spectrometer.


Copyright ©2019, Edinburgh Instruments Ltd. All rights reserved.

Sponsored by



### More Application Notes from this Sponsor

- Electroluminescence and Photoluminescence Spectroscopy of a Phosphorescent Organic Light Emitting Diode (PhOLED)
- Photophysical Characterisation of Perovskite Quantum Dots
- Charge Carrier Recombination Dynamics of Semiconductor Photocatalyst



Visit Photonics Media to download other white papers and learn more about the latest developments in lasers, imaging, optics, biophotonics, machine vision, spectroscopy, microscopy, photovoltaics and more.

[www.photonics.com/WhitePapers.aspx](http://www.photonics.com/WhitePapers.aspx)

We respect your time and privacy. You are receiving this email because you are a Photonics Media subscriber, and/or a member of our website, Photonics.com. You may use the links below to manage your subscriptions or contact us.

Questions: [info@photonics.com](mailto:info@photonics.com)

[Unsubscribe](#) | [Subscribe](#) | [Subscriptions](#) | [Privacy Policy](#) | [Terms and Conditions of Use](#)

Photonics Media, 100 West St., PO Box 4949, Pittsfield, MA 01202-4949

© 1996 - 2019 Laurin Publishing. All rights reserved. Photonics.com is Registered with the U.S. Patent & Trademark Office. Reproduction in whole or in part without permission is prohibited.