

PHOTONICS spectra

WHITE PAPERS & APPLICATION NOTES

APPLICATION NOTE

Temperature Dependent Triplet States of Benzophenone: Spectral and Lifetime Measurements Utilising Transient Absorption Spectroscopy



APL_PNA_Media_1020

Introduction

Triplet excited states in molecules and materials play an important role in many practical applications, from phosphorescent pigments to photovoltaic dyes, and even photovoltaic solar cells. New materials that exhibit triplet state photochemistry are constantly under development, and understanding the lifetime and energy transfer processes of their triplet states is a fundamental step in their design and application. In this application note, we demonstrate how transient absorption spectroscopy and temperature-dependent measurements can be employed to investigate and understand the nature of photoexcited triplet states of a molecule.

Triplet states are electronic states with two unpaired electrons with parallel spins in different molecular orbitals. Closed shell molecules have a single state (opposite spin pair) in their ground state, with their triplet being at higher energy. Figure 1 shows a schematic diagram representing the different energy levels for a diatomic molecular dication and the mechanisms for triplet state generation and deactivation. The ground state S_0 has two electrons occupying its bonding orbital (σ). Excitation of one electron into the antibonding (σ^*) orbital takes the molecule into its first excited singlet state S_1 . Transitions between states of the same spin multiplicity (singlet-singlet or triplet-triplet) are allowed by selection rules, but singlet-triplet transitions are spin-forbidden. In some cases, the molecule can undergo intersystem crossing to a lower lying triplet (T_1). This is favoured by molecules containing heavy atoms and paramagnetic species in solution.

Mixtures of a triplet excited state can proceed via several routes. The molecule may decay to its ground singlet state emitting light, in a process known as phosphorescence. Intermolecular triplet-triplet annihilation without the assistance of catalysis may also take place. A further possibility is the reaction of the molecule with its oxidant environment inducing relaxation to S_0 or quenching, as shown in Figure 1 where oxygen is converted to its triplet state by a molecule in the T_1 state.

Another process, not depicted in Figure 1, is triplet-triplet annihilation (TTA), where two triplet states react with each other to give a single singlet excited state, $T_1 + T_1 \rightarrow S_1 + S_0$. TTA has many possible outcomes: S_1 can decay via fluorescence or internal conversion, or S_0 and S_1 can form a triplet excimer which then dissociates.

The different pathways discussed above make triplet photochemistry a complex subject to study. However, modern optical spectroscopy techniques exist at investigating these excited state processes. Photoluminescence (PL) spectroscopy is a commonly used technique to probe triplet excited state lifetimes, but it is only effective when there is significant phosphorescence emission from the $T_1 \rightarrow S_0$ transition. At more temperature non-radiative losses often outweigh the radiative phosphorescence decay pathway, so measurements at low temperature or using other techniques are required. An often used alternative is nanosecond transient absorption (ns-TA) spectroscopy, which probes the triplet ($T_1 \rightarrow T_2$) absorption as a function of time. By directly probing the concentration of the triplet as opposed to its luminescence, it is possible to study other decay pathways such as quenching and TTA processes. In addition, measurements in a cuvette may be performed to understand the temperature-dependent dynamics of the system.

In this application note we study the photo-generated triplet state of benzophenone by ns-TA and photoluminescence spectroscopy at different temperatures. Benzophenone (Figure 2) is an efficient triplet sensitizer thanks to its high (~100%) intersystem crossing yield. Its S_1 state is generated by promoting an electron from a non-bonding orbital to the σ^* orbital from the carbonyl group, hence the (π, π^*) label in Figure 2. The higher excited state S_2 arises from the excitation of σ to σ^* (σ, σ^*) state. Intersystem crossing from S_2 to T_1 is very efficient because transitions between (π, π^*) and (σ, σ^*) states are highly favoured. This results in a large concentration of triplet states which can evolve via phosphorescence, non-radiative relaxation or triplet-triplet annihilation (TTA) events.




Figure 1. Schematic diagram of excitation of benzophenone to the S_1 and S_2 states, the subsequent intersystem crossing (ISC) to the T_1 state, and the subsequent decay pathways.

Copyright © 2020 Edinburgh Instruments Ltd. All rights reserved.

Temperature Dependent Triplet States of Benzophenone

Triplet excited states in molecules and materials play an important role in many practical applications, from phosphorescent pigments to photodynamic therapy, and even photovoltaic solar cells. New materials that exhibit triplet state photochemistry are constantly under development, and understanding the lifetime and energy transfer processes of their triplet states is a fundamental step in their design and optimisation. In this application note, we demonstrate how transient absorption spectroscopy and temperature-dependent measurements can be employed to investigate and understand the nature of photoexcited triplet states of a molecule.

[DOWNLOAD WHITE PAPER](#)



EDINBURGH
INSTRUMENTS

More Application Notes from This Sponsor

- [Electroluminescence and Photoluminescence Spectroscopy of a Phosphorescent Organic Light Emitting Diode \(PhOLED\)](#)
- [Photophysical Characterisation of Perovskite Quantum Dots](#)
- [Charge Carrier Recombination Dynamics of Semiconductor Photocatalyst](#)

Visit [Photonics Media](#) to download other white papers and learn more about the latest developments in lasers, imaging, optics, biophotonics, machine vision, spectroscopy, microscopy, photovoltaics and more.

www.photonics.com/WhitePapers.aspx

We respect your time and privacy. You are receiving this email because you are a Photonics Spectra magazine subscriber. You may use the links below to manage your subscriptions or contact us.

Questions: info@photonics.com

[Unsubscribe](#) | [Subscribe](#) | [Subscriptions](#) | [Privacy Policy](#) | [Terms and Conditions of Use](#)

Photonics Media, 100 West St., PO Box 4949, Pittsfield, MA 01202-4949

© 1996 - 2020 Laurin Publishing. All rights reserved. Photonics.com is Registered with the U.S. Patent & Trademark Office.

Reproduction in whole or in part without permission is prohibited.



Laurin Publishing

PHOTONICS MEDIA