

# PHOTONICS spectra®

## WHITE PAPERS & APPLICATION NOTES

**APPLICATION NOTE**

**Comparing Hole Extraction Efficiencies in Perovskite Solar Cells using Photoluminescence Quantum Yield**



AN\_P16, Stuart Thomson

**Introduction**

Efficient extraction of charge carriers is vital for the creation of high-efficiency solar cells. To improve charge extraction, electron and hole extraction layers are routinely incorporated into the solar cell stack. To optimise cell performance, the extraction efficiency of different extraction layers must be investigated and compared. Photoluminescence (PL) spectroscopy is an ideal technique for this comparison as PL is proportional to the number of charge carriers in the perovskite.

For PL to occur, the photo-generated electrons and holes in the perovskite layer must recombine. The transfer of electrons or holes into adjacent layers will inhibit radiative electron-hole recombination within the perovskite layer and decrease the PL response (Figure 1). By monitoring the intensity of the PL, the efficiency the electrons or holes are extracted out of the perovskite layer can be studied.




Figure 1: Intensity between incident laser source, photoluminescence and emitted charge extraction (into adjacent layers).

One promising material under investigation as a hole extraction layer in perovskite solar cells are vertically aligned carbon nanotubes (VACNTs). In this application note, the extraction efficiency of a VACNT hole extraction layer is investigated using PL spectroscopy with an Edinburgh Instruments FL51000 Photoluminescence Spectrometer.

**Experimental Setup**

A layer of mixed halide  $\text{CH}_3\text{Pb}(\text{I}_{0.8}\text{MA}_{0.2})\text{Br}_2$  perovskite was spin-coated atop a Glass/ITO and a Glass/ITO/VACNT substrate. Full details of the sample fabrication can be found in the accompanying publication. The PL of the two perovskite samples was characterised using the Edinburgh Instruments FL51000 photoluminescence spectrometer. The FL51000 was equipped with a 650 W Xenon lamp, double excitation and double emission monochromators and a PMT-980 detector.

For the relative PL intensity comparison, the perovskite cells held using the NAB3 Front Face Sample Holder and for determination of the PL quantum yield, they were placed in the NAB01 Integrating Sphere.



Figure 2: FL51000 Photoluminescence Spectrometer.

**Relative PL Intensity Comparison**

The simplest method to investigate the impact of the VACNT layer on the hole extraction efficiency is by measuring the PL spectra of the Glass/ITO/Perovskite and Glass/ITO/VACNT/Perovskite and comparing their integrated PL intensities. If the addition of the VACNT layer enhances hole extraction the PL intensity from the VACNT layer sample will be lower than the sample without the layer (Figure 3).



Figure 3: Relative PL intensity comparison.

For the relative PL intensity comparison, the perovskite samples were held in turn using the NAB03 Front Face Sample Holder of the FL51000 and their PL spectra measured under identical conditions (integration time, sample position and angle, excitation wavelength, and excitation and emission bandwidth) and are shown in Figure 4. It can be seen that the perovskite sample containing the VACNT layer has a significantly reduced PL contribution. Integrating the PL spectra using the Fluascan® software of the FL51000 reveals that the

## Comparing Hole Extraction Efficiencies in Perovskite Solar Cells using PLQY

Efficient extraction of charge carriers is vital for the creation of high-efficiency solar cells. Photoluminescence is proportional to the number of charge carriers in the perovskite and sensitive to charge transfer into adjacent layers. This makes PL based techniques invaluable for investigating the performance of new extraction layers. In this application note, the extraction efficiency of a VACNT hole extraction layer is investigated using PL spectroscopy with a Photoluminescence Spectrometer.

[DOWNLOAD WHITE PAPER](#)

Copyright ©2021, Edinburgh Instruments Ltd. All rights reserved.



## More White Papers from This Sponsor

- X-ray Excited Luminescence Lifetime of Scintillators
- Time-resolved Spectroscopy of Phosphorescent Oxygen Sensors in a Relevant in-vitro Environment for Biomedical Applications
- Raman Spectroscopy as a Tool for Studying Polymer Phase Transitions

Visit [Photonics Media](#) to download other white papers and learn more about the latest developments in lasers, imaging, optics, biophotonics, machine vision, spectroscopy, microscopy, photovoltaics and more.

[www.photonics.com/WhitePapers.aspx](http://www.photonics.com/WhitePapers.aspx)

We respect your time and privacy. You are receiving this email because you are a Photonics Spectra magazine subscriber. You may use the links below to manage your subscriptions or contact us.

Questions: [info@photonics.com](mailto:info@photonics.com)

[Unsubscribe](#) | [Subscribe](#) | [Subscriptions](#) | [Privacy Policy](#) | [Terms and Conditions of Use](#)

Photonics Media, 100 West St., PO Box 4949, Pittsfield, MA 01202-4949

© 1996 - 2022 Laurin Publishing. All rights reserved. Photonics.com is Registered with the U.S. Patent & Trademark Office. Reproduction in whole or in part without permission is prohibited.



LAURIN PUBLISHING

PHOTONICS MEDIA