

PHOTONICS spectra®

WHITE PAPERS & APPLICATION NOTES

TECHNICAL NOTE

Fluorescence Microscopy with the F55 Spectrofluorometer
TN_P55_Marco Fariña



Introduction

Fluorescence microscopy is an imaging technique that uses fluorescence, either by intrinsic emission or artificially added fluorophores, to provide contrast to microscope images. In biomedical imaging it offers a high degree of specificity and selectivity thanks to the use of multiple fluorescence labeling where different fluorophores can be used to simultaneously identify specific target molecules and cellular structures. Whereas, in materials science, fluorescence (photo-luminescence) microscopy can be used to observe photonic nanoparticles and image the surface of semiconductors for defects.




Figure 1 Composite (Widefield) fluorescence image of bovine pulmonary artery endothelial (BPAE) cells acquired using the F55 Spectrofluorometer with microscope add-on.

Fluorescence microscopy analysis detailed cellular images like that of bovine pulmonary artery endothelial (BPAE) cells shown in Figure 1. BPAE cells contain angiotensin-converting enzyme (ACE), which is a vital component in the regulation of blood pressure by constricting and dilating blood vessels. In the BPAE sample, the different cell components have been labelled with various colour fluorescent dyes in order to observe the structure of the cells. The technical note shows how a fluorescence image like that of BPAE above, along with emission spectra and lifetime decays, can be acquired using an Edinburgh Instruments F55 Spectrofluorometer with the microscope add-on.

Experimental Setup

The sample investigated was FluorGels™ Prepared Slide #1 from Invitrogen™ which contains bovine pulmonary artery endothelial (BPAE) cells. The BPAE cells are stained with a combination of fluorescent dyes, each targeting a specific structure in the cell. MuTrack™ Real DMFPlus is used to give red emission to mitochondria, Alexa Fluor™ 488 phalloidin to give green emission to the filamentous actin network, and blue emitting DAPI to label the nuclei.




Figure 2 F55 Spectrofluorometer coupled to a fluorescence microscope.

All measurements were acquired using an Edinburgh Instruments F55 Spectrofluorometer coupled to a Nikon N4U Upright Fluorescence Microscope. The F55 was fitted with the SC50 Liquid Light Guide Launcher and the excitation and emission light were coupled to and from the microscope using liquid light guides. The F55 was equipped with a 150W Xenon lamp for steady-state excitation, Time-Correlated Single Photon Counting (TCSPC) Module electronics, and a PMT-980 detector. The microscope was equipped with a CMOS fluorescence imaging camera for widefield fluorescence imaging, an age fluorescence excitation light guide coupler, an emission light guide coupler, and a laser excitation coupler with an EP-275 pulsed diode laser.

Copyright ©2022, Edinburgh Instruments Ltd. All rights reserved.

Spectrofluorometer Add-on for Widefield Fluorescence Microscopy

Fluorescence microscopy is an imaging technique that uses fluorescence, either by intrinsic emission or artificially added fluorophores, to provide contrast to microscope images. This technical note shows how fluorescence images along with emission spectra and lifetime decays, can be acquired using a Spectrofluorometer with the microscope add-on.

[DOWNLOAD WHITE PAPER](#)



More Application Notes from This Sponsor

- [Comparing Hole Extraction Efficiencies in Perovskite Solar Cells using PLQY](#)
- [X-ray Excited Luminescence Lifetime of Scintillators](#)
- [Time-resolved Spectroscopy of Phosphorescent Oxygen Sensors in a Relevant in-vitro Environment for Biomedical Applications](#)

Visit [Photonics Media](#) to download other white papers and learn more about the latest developments in lasers, imaging, optics, biophotonics, machine vision, spectroscopy, microscopy, photovoltaics and more.

www.photonics.com/WhitePapers.aspx

We respect your time and privacy. You are receiving this email because you are a Photonics Spectra magazine subscriber. You may use the links below to manage your subscriptions or contact us.

Questions: info@photonics.com

[Unsubscribe](#) | [Subscribe](#) | [Subscriptions](#) | [Privacy Policy](#) | [Terms and Conditions of Use](#)

Photonics Media, 100 West St., PO Box 4949, Pittsfield, MA 01202-4949

© 1996 - 2022 Laurin Publishing. All rights reserved. Photonics.com is Registered with the U.S. Patent & Trademark Office.

Reproduction in whole or in part without permission is prohibited.



Laurin Publishing

PHOTONICS MEDIA