

PHOTONICS spectra®

WHITE PAPERS & APPLICATION NOTES

Evaluation Warning: This version can be used only for evaluation purposes.

APPLICATION NOTE

Relative Quantum Yield of 2-Aminopyridine

AN_P18, Laura Nather, Daniela Gokarny, Georgios Avramitis



Introduction

Quantum yield is a fundamental photophysical parameter that describes a sample's fluorescence efficiency and is defined as the ratio of the number of photons emitted to the number of photons absorbed by a sample. Accurate and reliable quantum yield measurements are important for a broad range of applications including displays, solar cells, bioimaging and drug development.

There are two optical methods for measuring the quantum yield: the absolute method and the relative method. In the absolute method, the quantum yield is measured directly using an integrating sphere, while in the relative method, the fluorescence intensity of the unknown sample is compared with the fluorescence intensity of a standard sample to calculate the quantum yield of the unknown.

In this application note, an Edinburgh Instruments F55 Spectrofluorometer (Figure 1) is used to measure the quantum yield of 2-Aminopyridine (2AMP) via the relative method. 2AMP in sulfuric acid (H₂SO₄) has been previously used as a quantum yield reference standard in the UV-visible range. The quantum yield of 2AMP was measured to be 0.25 in 1961¹ and 0.16 in 1965.² These literature quantum yield reference values are now decades old, and here we present a re-investigation and evaluation of the quantum yield of 2AMP in 1M H₂SO₄ using quinine bisulfate (QBS) in 1M H₂SO₄ as the reference standard with a modern spectrofluorometer.



Figure 1: Edinburgh Instruments F55 Spectrofluorometer

Methodology

The relative quantum yield of 2AMP can be calculated through the following formula:³

$$\Phi_2 = \Phi_1 \left(\frac{I_2}{I_1} \right) \left(\frac{1 - 10^{-A_1}}{1 - 10^{-A_2}} \right) \left(\frac{n_1}{n_2} \right)^2 \quad (Eq. 1)$$

where the subscripts S and R denote samples of interest (2AMP) and reference standard (QBS), respectively. Φ is the quantum yield, I is the integrated fluorescence intensity, and A is the

Copyright ©2022, Edinburgh Instruments Ltd. All rights reserved.

absorbance at the excitation wavelength, n is the refractive index of the solvent, used for sample and reference solutions at the mean emission wavelength. In this application note, the same solvent (1M H₂SO₄) was used for both 2AMP and QBS and this term cancels.

To increase the accuracy and precision of the calculated quantum yield value, it is best practice to prepare and measure several solutions of the sample and reference with different concentrations. By plotting $\log(\Phi_2)$ for 2AMP and QBS, the gradients (G_{2AMP} and G_{QBS}) can be used to calculate the quantum yield (Eq.2). This approach prevents potential systematic errors, such as dye aggregation, which would appear as a deviation from the straight line at higher concentrations.

$$\Phi_2 = \Phi_1 \left(\frac{G_2}{G_1} \right) \left(\frac{1 - 10^{-A_1}}{1 - 10^{-A_2}} \right)^2 \quad (Eq. 2)$$

$$\Phi_2 = \Phi_1 \frac{G_{2AMP}}{G_{QBS}}$$

Five solutions of 2AMP in 1M H₂SO₄ and five solutions of QBS in 1M H₂SO₄ at different concentrations were prepared. Absorption and fluorescence spectra were measured using an F55 Spectrofluorometer equipped with a 150 W Xenon lamp, a PMT-PBS detector and a SC-05 cuvette holder.

Absorption & Emission Spectra of 2AMP & QBS

Firstly, the absorbance values of the five 2AMP and QBS solutions were determined by measuring the absorption spectra using the F55's built-in transmission detector. The absorbance values of the solutions were kept below 0.1 at the excitation wavelength (210 nm) to minimize the probability of inner filter effects and ranged between 0.009 and 0.095. The normalized absorption spectra of 2AMP and QBS are shown in Figure 2a.

Next, the fluorescence spectra of five 2AMP and QBS solutions were acquired. The intensity of the fluorescence detected by a spectrofluorometer is dependent on the excitation wavelength, excitation and emission bandwidths, and integration time. By keeping these parameters identical the integrated fluorescence intensities, I_2 and I_1 , of 2AMP and QBS can be meaningfully compared. The experimental parameters were $\lambda_{exc} = 210$ nm with excitation and emission bandwidths set at 5 nm and 0.5 nm, respectively, 1 nm step size, and 0.5 s integration time. Figure 2b shows the normalized fluorescence spectra of 2AMP and QBS.

Relative Quantum Yield of 2-Aminopyridine

Quantum yield is a photophysical parameter describing a sample's fluorescence efficiency. It can be measured optically via the absolute method and relative method. 2-Aminopyridine in sulfuric acid has been previously used as a QY reference standard in the UV-vis range but the 2 AMP QY values are now decades old. Edinburgh Instruments presents a reinvestigation and reevaluation of the quantum yield of 2 AMP in 1 M H₂SO₄, using quinine bisulphate in 1 M H₂SO₄ as the reference standard with a Spectrofluorometer.

[DOWNLOAD WHITE PAPER](#)



More White Papers from This Sponsor

- Spectrofluorometer Add-on for Widefield Fluorescence Microscopy
- Comparing Hole Extraction Efficiencies in Perovskite Solar Cells using PLQY
- X-ray Excited Luminescence Lifetime of Scintillators

Visit [Photonics Media](#) to download other white papers and learn more about the latest developments in lasers, imaging, optics, biophotonics, machine vision, spectroscopy, microscopy, photovoltaics and more.

www.photonics.com/WhitePapers.aspx

We respect your time and privacy. You are receiving this email because you are a Photonics Spectra magazine subscriber. You may use the links below to manage your subscriptions or contact us.

Questions: info@photonics.com

[Unsubscribe](#) | [Subscribe](#) | [Subscriptions](#) | [Privacy Policy](#) | [Terms and Conditions of Use](#)

Photonics Media, 100 West St., PO Box 4949, Pittsfield, MA 01202-4949

© 1996 - 2022 Laurin Publishing. All rights reserved. Photonics.com is Registered with the U.S. Patent & Trademark Office. Reproduction in whole or in part without permission is prohibited.



Laurin Publishing

PHOTONICS MEDIA