

PHOTONICS spectra®

WHITE PAPERS & APPLICATION NOTES



Precision Motion And How to Achieve It

PRECISION MOTION REQUIRES FASTIDIOUS ATTENTION TO INSTRUMENTATION, OTHER PASSIVE CONTROLS, AND EXPERIMENT DESIGN

By Jerica Cox Vog, Ph.D., Mad City Labs, Inc.

Achieving precision motion at the micron level is an impressive feat, a carefully choreographed mix of experiment design, instrumentation, and a keen eye for detail. It starts to reason that progressing to nanoscale-level precision would simply be a matter of refining a few steps. However, controlling motion at the nanoscale is a completely different dance.

Control an experiment in optical microscopy. A researcher may think, justifying that their micron-level experiment setup and instrumentation allow for precise control. Simply adding a nanopositioner will allow them to go the last mile in terms of motion control to take thinner slices of their samples and to locate more minute areas of interest within the sample. But researchers may not realize that sample image noise are not the only thing magnified; minor noise at the micron scale can be devastating at the nanometer scale.

Suddenly, things are not as frictionless as expected and everything seems to be moving. In fact, everything probably has been moving the whole time; the researcher just lacked the technical capability to observe that movement. In this respect, progressing from the micron scale to the nanoscale is like wearing a new pair of glasses.

What Factors Affect Precision Motion at Nanometer Scale?

A common misconception surrounding nanoscale motion is that precision depends solely on instrumentation — specifically, securing equipment with a technical specification of nanometer resolution. But any experiment's methodology must also account for thermal expansion, force distortion, misalignment, and other passive factors that affect both instrumentation and the experiment. The instrumentation a researcher chooses is part of passive control, while passive factors include acoustic effects and thermal centers of a workspace.

Acoustic effects are composed of any source of vibration, from people scanning lab doors to a distant meeting in the room above the lab. Thermal control dictates that researchers maintain optimal temperatures across the lab, down to the immediate area surrounding instrumentation.

Admittedly, it is easy to say, "you need to control your temperature and the acoustics." But a disparity usually exists between

the control researchers can exercise over their laboratory and the control they can exert over an entire building. Some researchers may be in a position to control items regarding their environment or workspace but, more often than not, they have been afforded a space and instructed to make a lab of it.

Clearly, researchers always should use the highest-quality instrumentation and data acquisition systems possible (e.g., passive controls). Beyond that, understanding what to look out for and recognizing their effects at the nanoscale allows experimenters to make the best of their respective situations. Researchers use almost university that they do not want to experiment near their spaces and that the room needs to be as quiet (i.e., vibration-free) as possible. But they may not understand the extent to which actions they perform as part of their workflow produce disruptive heat or vibration.

Embracing a Nano Mindset

For most Mad City Labs customers that have trouble achieving precision motion at the nanoscale, technical support does not comprise giving them new information. It is more about properly contextualizing information they already have, explaining how it affects their nanoscale experiment and how those effects can be mitigated. The attempt to bridge the gap between researchers' stated goals and achieving those goals in a practical sense, acknowledging that many of our customers' experiments exist at the cutting edge of what is technically feasible.

These researchers are experts in their disciplines and their experiments, with a deep understanding of what they are trying to probe. We are experts in the application of instrumentation setup to achieve those research goals. Collaborative, in-depth discussion (deconstructing every part of the experiment, from the methodology and equipment to passive factors), often are required to identify issues and collaborate on a solution. Consider the following examples:

Instrumentation — A researcher may be using a motorized nanopositioning stage with a nanopositioner set atop it. The nanopositioning stage is set to execute a same step that constantly measures its position and then retracts accordingly. While that may sound like a great idea, constantly sensing at a microscopic level has a huge impact on the microscope, potentially negating all the benefits of using a nanopositioner. To ensure nanopositioner effectiveness, the nanopositioner must be assigned a location and hold it within noise.

Precision Motion and How to Achieve It

Nanoscale motion requires fastidious attention to instrumentation, passive controls, and experiment design. Most users have experience in motion control within the sphere of the microscale. In this paper we discuss what needs to be considered in the jump to the nanoscale, and how detailed and deliberate choices at the design phase can lead to successful outcomes.

[DOWNLOAD WHITE PAPER](#)



More White Papers from This Sponsor

- Ingenuity + Instrumentation: Creating a Novel 4D Microscopy Method
- Nanopositioning! What is it good for?
- Resonant probe AFM: Uses and Advantages

Visit [Photonics Media](#) to download other white papers and learn more about the latest developments in lasers, imaging, optics, biophotonics, machine vision, spectroscopy, microscopy, photovoltaics and more.

www.photonics.com/WhitePapers.aspx

We respect your time and privacy. You are receiving this email because you are a Photonics Spectra magazine subscriber. You may use the links below to manage your subscriptions or contact us.

Questions: info@photonics.com

[Unsubscribe](#) | [Subscribe](#) | [Subscriptions](#) | [Privacy Policy](#) | [Terms and Conditions of Use](#)

Photonics Media, 100 West St., PO Box 4949, Pittsfield, MA 01202-4949

© 1996 - 2022 Laurin Publishing. All rights reserved. Photonics.com is Registered with the U.S. Patent & Trademark Office. Reproduction in whole or in part without permission is prohibited.



LAURIN PUBLISHING

PHOTONICS MEDIA