

PHOTONICS spectra®

WHITE PAPERS & APPLICATION NOTES

APPLICATION NOTE

Monitoring Pesticides on Apple Skin Using SERS

AN_316_Angela Flack



Introduction

Protecting food crops is vital to the world's food chain, and pesticides are a critical tool to kill, repel, and control pests. Two million tonnes of pesticides are used annually worldwide, with this value predicted to increase over time. However, pesticides can cause severe environmental issues to aquatic systems, wildlife, air, and soil. Additionally, the reported effects of pesticides on human health range from irritation, such as skin irritation and headaches, to chronic effects, such as asthma and cancer. There are also concerns about the consequences of consuming pesticides from the soil, but repetitive, close to our diet.¹

One method to reduce the dangers of pesticides is limiting the allowed concentration used on crops. This is controlled by analyzing the crop and determining the level of pesticide present. Apples are one of the most pesticide-treated fruits used as such producers have a list of regulations they must abide by for sale to the public.




Figure 1 Chemical structure of phosmet

Phosmet is an organophosphate insecticide used on apple trees to control codling moths, the 2018 Code of Federal Regulations set the pesticide residue tolerance on apples of phosmet at 10 ppm. High-performance liquid chromatography, mass spectrometry, and gas chromatography are currently the most used techniques for testing down to low concentrations. However, these are time-consuming and costly, and Raman spectroscopy offers several advantages such as rapid and non-destructive, fingerprint-like identification with little to no sample preparation.

Surface-enhanced Raman scattering (SERS) is an enhancement technique used with Raman spectroscopy to provide lower limits of detection. SERS offers a signal enhancement of up to 10⁷-10¹¹ and also advantageously quenches the fluorescence of analytes. Commercial SERS substrates are readily available, making SERS an accessible technique for low-concentration detection. In this application note, residual phosmet insecticide on apple skin is detected using SERS.

Materials and Methods

An RNS Raman Microscope equipped with a 785 nm laser and a 600 g/mm grating was used to analyze the phosmet samples. Gold SERS substrates from Horiba were used to provide the enhancement effect. Since gold nanoparticles typically absorb in the region of 520 nm - 580 nm the preferred excitation wavelength for gold SERS substrates is standardly 785 nm.



Figure 2 Experiment setup for apple analysis

Aqueous solutions of 10 ppm and 1 ppm phosmet were prepared and 6 µL of each phosmet solution was pipetted into the well of the SERS substrate, left to dry, and analyzed. To create a contaminated apple sample, a small volume of 10 ppm phosmet solution was dropped onto a cleaned, cut-out section of apple skin and left to dry. Once dry a wet cotton bud was used to smudge the surface of the apple skin and placed into a small volume of water. 6 µL of the sample was pipetted onto the SERS substrate for analysis.

The SERS substrates were analyzed by placing underneath the Raman microscope and exciting the sample with a 785 nm laser using an exposure time of 5 seconds. As no spatial resolution was required, the pinhole diameter was set at 2 mm to maximize the Raman signal.

Results and Discussion

Figure 3 shows the Raman spectra of the 4 samples analysed. The 10 ppm solution was analysed in a cuvette meaning there was no SERS enhancement; this measurement was performed to show the standard Raman response from phosmet at a low concentration. Without an enhancement from the SERS substrate, the 10 ppm solution of phosmet has very little Raman signal; however, when deposited on the gold SERS substrate a strong Raman spectrum is observed. When further diluted to 1 ppm a pronounced Raman signal is still observed, exceeding the low detection limit achievable with SERS. The spectrum shown in pink in Figure 3 is the washed apple skin spiked with 10 ppm phosmet. The spectrum is noisier than the 10 ppm solution of phosmet directly on the SERS substrate however well-defined peaks are still observed, peak assignments are detailed in Figure 4.

Pesticide Detection on Apple Skin Using SERS

Protecting food crops is vital to the world's food chain, and pesticides are a critical tool to kill, repel, and control pests. One method to reduce the dangers of pesticides is limiting the allowed concentration used on crops. This is controlled by analyzing the crop and determining the level of pesticide present. Apples are one of the most pesticide-treated fruits and as such producers have a list of regulations they must abide by for sale to the public. In this application note, find out how residual phosmet insecticide on apple skin is detected using Surface Enhanced Raman Scattering (SERS).

[DOWNLOAD WHITE PAPER](#)



More White Papers from This Sponsor

- [Multiphoton Imaging of Mouse Intestine](#)
- [Whisky Analysis by Raman Spectroscopy](#)
- [Measuring Picosecond Fluorescence Lifetimes Using the FLS1000 Equipped with a Hybrid Photodetector](#)

Visit [Photonics Media](#) to download other white papers and learn more about the latest developments in lasers, imaging, optics, biophotonics, machine vision, spectroscopy, microscopy, photovoltaics and more.

www.photonics.com/WhitePapers.aspx

We respect your time and privacy. You are receiving this email because you are a Photonics Spectra magazine subscriber. You may use the links below to manage your subscriptions or contact us.

Questions: info@photonics.com

[Unsubscribe](#) | [Subscribe](#) | [Subscriptions](#) | [Privacy Policy](#) | [Terms and Conditions of Use](#)

Photonics Media, 100 West St., PO Box 4949, Pittsfield, MA 01202-4949

© 1996 - 2023 Laurin Publishing. All rights reserved. Photonics.com is Registered with the U.S. Patent & Trademark Office. Reproduction in whole or in part without permission is prohibited.



Laurin Publishing

PHOTONICS MEDIA