

PHOTONICS spectra®

WHITE PAPERS & APPLICATION NOTES

APPLICATION NOTE
Plastics Identification Using ATR-FTIR Spectroscopy

AN_002, George Avontzakis



Introduction

Plastics are among the most used products worldwide, with global plastics production in 2020 estimated to be 367 million metric tons.¹ Plastics are incredibly versatile due to their low-cost production, durability, lightweight, and design freedom. As such, they can be found in almost every sector. Their wide range of applications spans from the food industry (for packaging) to the electrical and electronics sector (for insulation) and the construction industry for roofing, cladding, etc.

However, the massive production of plastics has led to a throw-away culture rendering plastic pollution one of the most ubiquitous environmental issues. One solution to prevent plastic pollution is improving waste management systems and recycling. As such, separating plastics based on their type is essential to ensure their suitable disposal and reuse. To facilitate plastics separation, the Society of the Plastics Industry introduced the Plastics Identification Code (PIC), which can usually be seen on the surface of plastics. Table 1 shows the PIC symbols with the acronyms, full names, and application examples. When a PIC identifier is absent, spectroscopic techniques, such as Fourier Transform Infrared Spectroscopy (FTIR), can quickly identify and sort plastics for recycling.

Table 1. PIC symbols for plastics with their acronyms, full names, and application examples.

Symbol	Acronym	Full Name	Application Examples	Recyclable
	PET	Polylethylene Terephthalate	Plastic bottles, soft-drink cans	✓
	HDPE	High-Density Polyethylene	Milk crates, water bottles, pipe, shower stalls	✓
	PPC	Polycarbonate	Optical lenses, safety glasses	✗
	LDPE	Low-Density Polyethylene	Plastic food wrap, household bottles, milk jugs	✓
	PE	Polypropylene	Plastic toys, chairs, crates, shower seats	✓
	PS	Polystyrene	CD cases, plastic cups	✗
	Other	Other plastics	Plastic containers, packaging materials, medical devices	✗

FTIR is a powerful qualitative and quantitative analysis tool. When a plastic absorbs IR radiation, the resulting signal is a spectrum that represents its molecular "fingerprint". Different plastic samples produce different fingerprints that enable

plastic identification. In this application note, an Edinburgh Instruments IR FTIR Spectrometer equipped with an Attenuated Total Reflection (ATR) accessory was used to identify different plastics.

Materials and Methods

Four plastic samples were analyzed: a biodegradable bag, a plant-based bag for electronics packaging, a soft-drink bottle, and a juice bottle. The four types were cut to roughly 2x2 cm² in size. The samples were analyzed using an Edinburgh Instruments IR FTIR Spectrometer equipped with an ATR accessory with a ZnSe crystal. The ATR accessory requires minimal to no sample preparation and enables the acquisition of reliable and high-quality spectra. To obtain a spectrum, the clamping arm of the ATR accessory applies pressure on the sample, ensuring consistent contact between the crystal and the sample. The plastics were identified from their IR spectra using the KnowItAll® spectral library.



Figure 1. Edinburgh Instruments IR FTIR Spectrometer.

Results and Discussion

Spectra of the four plastics were acquired with a resolution of 4 cm⁻¹, averaging 5 scans with a total acquisition time of around 15 seconds for each spectrum. The spectral library identified the biodegradable bag sample (Figure 2) as high-density polyethylene (HDPE) with a 90% match. HDPE is a "linear" form of PE whose polymer chains are absent of methyl (CH₃) side chains, containing only CH₂ groups. In such a way, the polymer chains can be packed closer together, resulting in a denser material. The band at 1465 cm⁻¹ is an HDPE's characteristic peak. The bands at 2915 cm⁻¹, 2848 cm⁻¹, and 719 cm⁻¹ correspond to the methylene (CH₂) group vibrations attributed to the asymmetric C-H stretch, the symmetric C-H stretch, and the rocking mode, respectively.¹¹ From Table 1, it is seen that the HDPE sample is a PIC 2 material and recyclable.

Plastics Identification Using ATR-FTIR Spectroscopy

The massive production of plastics has led to a throw-away culture rendering plastic pollution one of the biggest environmental issues. One solution to prevent plastic pollution is improving waste management systems and recycling. To facilitate plastics separation, the Plastics Identification Code (PIC) was introduced. When a PIC identifier is absent, spectroscopic techniques, such as Fourier Transform Infrared Spectroscopy (FTIR), can quickly identify and sort plastics for recycling.

[DOWNLOAD WHITE PAPER](#)



More White Papers from This Sponsor

- [ATR-FTIR of Blood Serum Using a Heated ATR Accessory](#)
- [UV-Vis Spectroscopy for Characterising the Optical Properties of Gold Nanoparticles](#)

Visit [Photonics Media](#) to download other white papers and learn more about the latest developments in lasers, imaging, optics, biophotonics, machine vision, spectroscopy, microscopy, photovoltaics and more.

www.photonics.com/WhitePapers.aspx

We respect your time and privacy. You are receiving this email because you are a Photonics Spectra magazine subscriber. You may use the links below to manage your subscriptions or contact us.

Questions: info@photonics.com

[Unsubscribe](#) | [Subscribe](#) | [Subscriptions](#) | [Privacy Policy](#) | [Terms and Conditions of Use](#)

Photonics Media, 100 West St., PO Box 4949, Pittsfield, MA 01202-4949

© 1996 - 2023 Laurin Publishing. All rights reserved. Photonics.com is Registered with the U.S. Patent & Trademark Office. Reproduction in whole or in part without permission is prohibited.