

BioPhotonics

Bringing Light to the Life Sciences®

WHITE PAPERS & APPLICATION NOTES

APPLICATION NOTE

Olive Oil Quality Assessment with UV-Vis Spectrophotometry

AN_UVO_Maria Tere



Introduction

Olive oil has been valued for its culinary and cultural significance for millennia. As one of the most widely consumed edible oils globally, quality standards are defined and monitored by regulatory bodies such as the International Olive Council (IOC). The IOC has defined grades for olive oil based on its quality and purity, with three primary classifications: extra virgin, virgin, and olive oil. Extra virgin olive oil is obtained solely by cold mechanical extraction and is prized for its superior taste, characterized by low acidity. Virgin olive oil, while slightly lower in quality than extra virgin, still maintains significant sensory and nutritional qualities. Olive oil, the most widely available grade, may undergo refining processes and often exhibits a milder flavor profile.

Despite various countermeasures, several cases of olive oil frauds are reported each year in Europe.¹ These cases commonly involve mixing extra virgin olive oil with lower-quality edible oil or adding colorants to the olive oil to pass as extra virgin. The integrity of olive oil can also be compromised by improper storage and processing. Reliable olive oil analysis methods are essential for ensuring adherence to regulatory standards and gaining consumer trust.

One of the most reliable techniques for assessing the quality of olive oil is UV-visible (UV-Vis) spectrophotometry. This analytical method offers a precise and efficient means of quantifying various compounds present in olive oil, including pigments, chemical compounds, and oxidation products, by measuring the absorbance of specific wavelengths of light. UV-Vis spectrophotometry only requires a benchtop instrument and does not require sample preservation; therefore, it is an ideal tool for the rapid identification of adulteration or degradation.

In this Application Note, an Edinburgh Instruments DSS Spectrophotometer is used for the quality assessment of olive oils, including extra virgin, virgin, and lower-grade edible olive oils, following the industry's standard UV-Vis analysis method. Additionally, a method to quantify olive oil adulteration using its visible absorption spectrum is evaluated.

Materials and methods

Commercial edible oils were purchased and used as provided (Figure 1, Table 1). Sample HEVOO was prepared by heating EVOO above its smoke point for a duration of 5 minutes.

Table 1: Table of sample names and their abbreviations.

Sample	Details
EVOO	Extra virgin olive oil
VVO	Virgin olive oil
PO	Olive pomace oil (blend of refined olive pomace oil and virgin olive oil)
HEVOO	Heated extra virgin olive oil (EVOO sample heated above smoke point for 5 minutes)
CO	Blend of oil: Commercial blend of 90% refined sunflower oil and 10% extra virgin olive oil
SO	Sunflower oil




Figure 1: Labeled sample containers and their abbreviations.

The quality of samples (EVOO, VVO, PO, and HEVOO) was assessed using the method established by the International Olive Council.² The method requires a spectrophotometer which can measure with 1 nm resolution between 220 nm and 960 nm. Linearity, wavelength and photometric accuracy in the UV range must be verified before the measurements. This study employed a DSS Spectrophotometer (Figure 2). Its performance was validated using halogen oxide and petroleum dicolorimetric standards. The samples were diluted in hexane to a concentration of 1% (w/v), transferred to quartz cuvettes and studied in the DSS Spectrophotometer after performing a baseline correction.



Figure 2: Edinburgh Instruments DSS UV-Vis Spectrophotometer.

Copyright © 2023, Edinburgh Instruments Ltd. All rights reserved.



More White Papers from This Sponsor

- Photoluminescence and Electroluminescence Confocal Imaging of an OLED
- Characterization of SERRS Nanoparticles Using UV-Vis and Raman Spectroscopy
- Plastics Identification Using ATR-FTIR Spectroscopy

Visit [Photonics Media](#) to download other white papers and learn more about the latest developments in lasers, imaging, optics, biophotonics, machine vision, spectroscopy, microscopy, photovoltaics and more.

www.photonics.com/WhitePapers.aspx

We respect your time and privacy. You are receiving this email because you are a Photonics Spectra magazine subscriber. You may use the links below to manage your subscriptions or contact us.

Questions: info@photonics.com

[Unsubscribe](#) | [Subscribe](#) | [Subscriptions](#) | [Privacy Policy](#) | [Terms and Conditions of Use](#)

Photonics Media, 100 West St., PO Box 4949, Pittsfield, MA 01202-4949

© 1996 - 2024 Laurin Publishing. All rights reserved. Photonics.com is Registered with the U.S. Patent & Trademark Office. Reproduction in whole or in part without permission is prohibited.



LAURIN PUBLISHING

PHOTONICS MEDIA