

# PHOTONICS spectra®

## WHITE PAPERS & APPLICATION NOTES



### Resonant Probe AFM: Uses and Advantages

RESONANT PROBE AFM: MCL RECOGNIZES UNIQUE CONSTRUCTION AND FUNCTIONALITY ENABLE THEM TO SERVE IN APPLICATIONS DEMANDING HIGH ACCURACY, ADAPTABILITY, AND EXCEPTIONALLY HIGH RESOLUTION

By James Mackay, Ph.D., Mad City Labs, Inc.

A resonant probe atomic force microscope (AFM) offers capabilities and benefits not found in commonly used optical deflection AFMs in the measurement of physical constants, physical properties, and surface morphology.

However, resonant probe AFMs' functionality and potential remain relatively unknown, even though resonant probe technology and its methodology are decades old. This article discusses resonant probe AFM construction, functionality, and typical applications, as well as introduces readers to MacCity Labs' custom AFM solutions.

#### Resonant Probe Construction And Functionality

The basis for resonant probe microscopy — and its key differentiator from optical deflection apparatus — is a piece of quartz crystal shaped like a tuning fork.

The technology was developed and leveraged by Sisko in the 1960s as an oscillator that would vibrate at a specified frequency, very accurately. The technology was applied by Sisko to timepieces, and later it was integrated into microprocessors to guide their "time of day" courses.

In the 1980s, IBM Research developed resonant probe microscopy using these quartz tuning forks at its Armonk innovation lab in Silicon Valley. Under this technique, the tuning fork — with its tines made of various materials applied to the end, depending on the application — is used as a force sensor.

The fundamental difference between the optical deflection technique and the resonant probe technique is that the resonant probe always is oscillating when operating. It is always in intermittent contact with the surface. Optical deflection systems were initially envisioned only as contact probes, wherein you drag the probe across a surface, absent any oscillation, and you can measure how it moves (e.g., measuring frictional forces on a surface). That said, intermittent contact probing was developed later for optical deflection AFMs.

Because of the extent to which resonant probe microscopy allows experimenters to modify and customize probe tips, it tends to attract individuals and teams who require customization and the flexibility to make their own tips.



At the advent of commercial AFM microscopy nearly 30 years ago, vendors used to produce standardized silicon-coated tips — essentially, small cantilevers wherein you drive a laser into the back of the tip, place that on a position sensor detector, and then measure the deflection of that tip. While the methodology has dominated the market since, the tuning fork has emerged as a suitable AFM microscopy tool.

For example, attaching a standard commercial optical deflection tip to the tuning fork allows it to be used for both optical deflection and resonant probe microscopy. A wire glued to the end of the probe could be used in one of two ways: as a conductive needle for direct imaging, and the wire tip could be used as a near-field probe; in, one could apply a bias voltage to the tip.

## Resonant probe AFM: Uses and Advantages

A resonant probe atomic force microscope (AFM) offers capabilities and benefits not found in commonly used optical deflection AFMs in the measurement of physical constants, physical properties, and surface morphology. This article discusses resonant probe AFM construction, functionality. Some of the applications particularly well suited to resonant probe AFM are optical antenna properties, quantum materials (NV center research), and terahertz microscopy/spectroscopy.

[DOWNLOAD WHITE PAPER](#)



## More White Papers from This Sponsor

- [Near-Field Scanning Optical Microscopes: Capabilities and Applications](#)
- [Atomic Force Microscopy Enhances the Nanoscopy Toolkit](#)
- [Understand Noise at the Sub-nanometer Level](#)

Visit [Photonics Media](#) to download other white papers and learn more about the latest developments in lasers, imaging, optics, biophotonics, machine vision, spectroscopy, microscopy, photovoltaics and more.

[www.photonics.com/WhitePapers.aspx](http://www.photonics.com/WhitePapers.aspx)

We respect your time and privacy. You are receiving this email because you are a Photonics Spectra magazine subscriber. You may use the links below to manage your subscriptions or contact us.

Questions: [info@photonics.com](mailto:info@photonics.com)

[Unsubscribe](#) | [Subscribe](#) | [Subscriptions](#) | [Privacy Policy](#) | [Terms and Conditions of Use](#)

Photonics Media, 100 West St., PO Box 4949, Pittsfield, MA 01202-4949

© 1996 - 2021 Laurin Publishing. All rights reserved. Photonics.com is Registered with the U.S. Patent & Trademark Office.

Reproduction in whole or in part without permission is prohibited.



LAURIN PUBLISHING

PHOTONICS MEDIA